



Centrifuge Tubes 15ml/50ml, ETO Sterile

Codes 563-819, 563-820



Centrifuge Tubes 15ml/50ml, ETO Sterile

Westlab's 15ml and 50ml Centrifuge Tubes are ideal for use in centrifugation, sample storage, and molecular biology.

Product Features

- Centrifuge Tubes, ETO Sterilised, Conical Bottom
- Medical Grade PP Tube & HDPE Leak-Proof Cap
- Non-Pyrogenic, Non-Cytotoxic & RNase/DNase Free
- Flat Top for Sample Identification Writing Area
- Resistant to Alcohols & Mild Organic Solvents
- Centrifugation up to 11,000 x g

Product Specifications

| Capacity | Sterile | Inner Bag | Outer Carton |
|----------|---------|------------|--------------|
| 50ml | Yes | 25 per bag | 500 |
| 15ml | Yes | 25 per bag | 1000 |

PRODUCT CODES

563-819
563-820

