



Steel Flammable Safety Cabinet

Code 555-5060, 555-5059, 555-5062, 555-5061,
555-5063, 555-5064, 555-5073



Steel Flammable Safety Cabinet

Product Features

- Double-wall 18-gauge welded steel construction with 1 1/2" air space
- Dual 2" capped vents with flame arrestors
- Full-height piano hinges open a full 180° for easy access
- Recessed handle
- 2" high leak proof door sill
- All shelves are adjustable on 2 1/2" centers
- Adjustable leveling feet
- Tough corrosive resistant finish
- Meets NFPA Code 30 and OSHA standard 1910.106 for storage for class I, II and III liquids

Capacity & Dimensions

555-5060	12 Gallon	Width 23 x Depth 18 x Height 35 Inches
555-5059	16 Gallon	Width 23 x Depth 18 x Height 44 Inches
555-5062	22 Gallon	Width 35 x Depth 22 x Height 35 inches
555-5061	30 Gallon	Width 43 x Depth 18 x Height 44 inches
555-5063	45 Gallon	Width 43 x Depth 18 x Height 65 Inches
555-5064	60 Gallon	Width 34 x Depth 34 x Height 65 Inches
555-5073	90 Gallon	Width 43 x Depth 34 x Height 66 Inches

*OSHA Code 29CFR 1910.106 indicates that no more than 60 gallons may be Class I, Class II or Class III flammable liquids.

PRODUCT CODE

555-5060
 555-5059
 555-5062
 555-5061
 555-5063
 555-5064
 555-5073

Canada

1 877 822 1455
 sales.ca@westlab.com
 www.westlab.com

Australia

1800 358 101
 sales.au@westlab.com
 www.westlab.com.au

New Zealand

+64 9 553 3677
 sales.nz@westlab.com
 www.westlab.co.nz

United States

646 568 5391
 sales.us@westlab.com
 www.westlab.com