

High RCF Centrifuge Tubes w/ Cap, 2 Sizes



These centrifuge tubes are engineered to meet the demands of the most critical scientific work. They are manufactured in a class 100,000 clean room, gamma irradiated to SAL 10⁻⁶ and are certified to be free from DNase, RNase, Human DNA, Pyrogens, and Endotoxins.

Specifications

- Ultra High Performance
- Leak Resistant
- Efficient automated production reduces cost and assures high quality
- Easy to use, clear tubes have black printed graduation marks with a bright white marking space

Item Code	Product Details
661-2311	Centrifuge Tubes w/cap 15 ml, Sterile, 17,000RCF, Case of 500, Bagged
661-2317	Centrifuge Tubes w/cap 50 ml, Sterile, 20,000RCF, Case of 500, Bagged