



Centrifuge Tubes 15ml/50ml, ETO Sterile

Codes 563-819, 563-820



Centrifuge Tubes 15ml/50ml, ETO Sterile

Westlab's 15ml and 50ml Centrifuge Tubes are ideal for use in centrifugation, sample storage, and molecular biology.

Product Features

- Centrifuge Tubes, ETO Sterilised, Conical Bottom
- Medical Grade PP Tube & HDPE Leak-Proof Cap
- Non-Pyrogenic, Non-Cytotoxic & RNase/DNase Free
- Flat Top for Sample Identification Writing Area
- Resistant to Alcohols & Mild Organic Solvents

Product Specifications

Capacity	Centrifugation	Inner Bag	Outer Carton
50ml	12,000 RPM (16,000 RCF)	25 per bag	500
15ml	8,000 RPM (11,000 RCF)	25 per bag	1000

PRODUCT CODES

563-819
563-820



Australia
1800 358 101
sales.au@westlab.com
www.westlab.com.au

Canada
1 877 822 1455
sales.ca@westlab.com
www.westlab.com

New Zealand
+64 9 553 3677
sales.nz@westlab.com
www.westlab.co.nz

United States
646 568 5391
sales.us@westlab.com
www.westlab.com