

# Westlab Laboratory Trolleys

Code 663-560 to 663-563



## Westlab Laboratory Trolleys

Westlab's range of Laboratory Trolleys have been designed for a multitude of applications. All the carts are extremely lightweight and have a safe working load of 250kg. They are supplied with precision bearing castors including 2 brakes.

Safe up to  
**250kg**

The construction of these units is from the highest quality, nylon fibre blended polymers which are resistant to the harshest working conditions.

We have switched the material to Recycled Polypropylene to reduce manufacturing costs which we have passed onto the customer.

### Product Features

- Heavy Duty Frame
- Manufactured from Recycled Polypropylene
- Removable Recessed Tops for Easy Cleaning
- Fitted into Solid Steel Blocks for Extra Strength
- 4 x 100mm Castors Supplied, 2 x Locked
- 250kg Load Rating

### PRODUCT CODES

663-560C  
663-561C  
663-562C  
663-563C

#### Australia

1800 358 101  
sales.au@westlab.com  
www.westlab.com.au

#### Canada

1 877 822 1455  
sales.ca@westlab.com  
www.westlab.com

#### New Zealand

+64 9 553 3677  
sales.nz@westlab.com  
www.westlab.co.nz

#### United States

646 568 5391  
sales.us@westlab.com  
www.westlab.com

## Product Specifications

Code	663-560C
Description	Small, 2 Shelf
Shelf Area	680 x 450mm
Packing Size	765 x 510 x 200mm H
Dimensions*	850 x 480 x 950mm H
Weight	11kg



Code	663-561C
Description	Small, 3 Shelf
Shelf Area	680 x 450mm
Packing Size	800 x 510 x 290mm H
Dimensions*	850 x 480 x 1020mm H
Weight	13.7kg



Code	663-562C
Description	Large, 2 Shelf
Shelf Area	920 x 520mm
Packing Size	1000 x 550 x 200mm H
Dimensions*	1100 x 520 x 950mm H
Weight	15kg



Code	663-563C
Description	Large 3 Shelf
Shelf Area	920 x 520mm
Packing Size	1000 x 600 x 290mm H
Dimensions*	1100 x 520 x 1020mm H
Weight	19kg



\* Dimensions include handles.

**Australia**  
1800 358 101  
sales.au@westlab.com  
www.westlab.com.au

**Canada**  
1 877 822 1455  
sales.ca@westlab.com  
www.westlab.com

**New Zealand**  
+64 9 553 3677  
sales.nz@westlab.com  
www.westlab.co.nz

**United States**  
646 568 5391  
sales.us@westlab.com  
www.westlab.com