


# Westlab Precision Balances

Codes 665-082C



## Westlab Precision Balances

Designed with high precision load cells and solid structures, Westlab's Precision Balances guarantee reliability and stability. The stainless-steel pans are corrosion resistant and easy to clean.

The integrated rechargeable battery allows the user to move from lab-to-lab or use the unit in an area with no access to a power source. In a teaching environment, science educators can charge the balances before a lesson to avoid the power cords being tangled. Additionally, the low or charged battery indicator displayed on the screen ensures that sample weighing is not interrupted.

With a backlit LCD screen, the balances feature a simple four key menu for easy navigation.

---

### Product Features

- LCD Display, Backlit
- High Precision Load Cell
- Integrated Rechargeable Battery
- Selectable Measuring Units
- Overload Display Indicator
- Piece Counting Function
- Height Adjustable Feet
- Intelligent Calibration
- Full Capacity Taring

#### Canada

1 877 822 1455

[sales.ca@westlab.com](mailto:sales.ca@westlab.com)

[www.westlab.com](http://www.westlab.com)

---



## Precision Balance 3000g x 0.01g

### Product Specifications

Display	LCD, backlit
Load Rating	3000g
Accuracy	0.01g
Tare Range	0 - 3000g
Repeatability	0.1g
Linearity	± 0.03g
Overload Protection	Yes
Feet	Height adjustable
Power Supply	220 – 240V or rechargeable battery
Pan Dimensions	200mm x 140mm
Unit Dimensions	220mm L x 220mm W x 70mm H

### PRODUCT CODES

665-082C

#### Canada

1 877 822 1455  
sales.ca@westlab.com  
www.westlab.com